BspQ I, HQ

G C T C T T C N
C G A G A A G N N N N

Code No. 1228A **Size:** **20,000 U**
Conc.: **20 U/μl**

Supplied Reagent:
10X BspQ I Buffer, HQ **10 ml**

Description:

BspQ I, HQ (high quality) is intended to be used for the preparation of active pharmaceutical ingredients in non-clinical trials and for the process development of drug manufacturing under GMP regulation. It can also be used for basic research of RNA therapeutics. The product is free from not only any human- or animal-derived materials but also β-lactam compounds in the final formulation.

This enzyme has an equivalent performance to BspQ I (Cat. #1227A/B), which recognizes asymmetric specific double-stranded DNA sequences and cleaves both strands outside the recognition sequence (Type IIS restriction enzyme). Therefore, in the preparation of linearized plasmids for *in vitro* transcription (IVT), introducing the recognition site of this enzyme immediately after the Poly(A) sequence of the target mRNA makes it possible to design mRNAs that do not have an extra sequence after the Poly(A) sequence.

Product Quality:

1. This product does not contain any human- or animal-derived materials in the final formulation.
2. This product does not contain any β-lactam compounds in the final formulation.

Storage Buffer:

20 mM	Tris-HCl, pH 7.5
500 mM	KCl
0.1 mM	EDTA
1 mM	DTT
0.1%	Tween 20
50%	Glycerol

10X BspQ I Buffer, HQ:

500 mM	Tris-HCl, pH 7.9
1 M	NaCl
100 mM	MgCl ₂
0.1%	Tween 20

Storage: -20°C

Source:

Escherichia coli carrying the plasmid containing the gene for BspQ I

Application:

Linearization of template plasmid used for mRNA synthesis

General Reaction Mixture and Incubation Time:

BspQ I, HQ*	1 μl
10X BspQ I Buffer, HQ	10 μl
Substrate DNA	≤ 10 μg
Nuclease-free water	up to 100 μl

↓

Incubate at 50°C for 1 hour.

* Add to the reaction lastly.

Reaction Temperature: 50°C

Note: The enzyme activity at 37°C decreases to about 30% of its activity at 50°C.

Precaution for Use:

Do not mix the enzyme vigorously.

Unit Definition:

One unit is defined as the amount of enzyme that cleaves 10 pmol of fluorescence-labeled probes in 1 hour at 50°C in a 50 μl reaction.

Heat Inactivation: 70°C for 20 min

Effect of DNA Methylation:

This enzyme is not sensitive to *dam*, *dcm*, and CpG methylation.

Star Activity:

Star activity may result from an excess incubation (e.g., > 4 hours).

Quality Control Data:

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

Related Products:

Template Vector (BspQ I) for T7 mRNA Synthesis (Cat. #6146)
In-Fusion® Snap Assembly Master Mix
(Cat. #638943/638944/638947 - 638949)

BspQ I (Cat. #1227A/B)
RNase-free Water (Cat. #9012)
Takara IVTpro™ T7 mRNA Synthesis Kit (Cat. #6144)
T7 RNA Polymerase, HQ (Cat. #2542A)
Pyrophosphatase (inorganic), HQ (Cat. #2451A)
PrimeCap™ T7 RNA Polymerase (low dsRNA) (Cat. #2560A)
T7 RNA Polymerase ver.2.0 (Cat. #2541A)
Pyrophosphatase (inorganic) (Cat. #2450A/B)
Recombinant RNase Inhibitor ver.2.0 (Cat. #2315A/B)
ATP (Cat. #4041)
GTP (Cat. #4042)
CTP (Cat. #4043)
UTP (Cat. #4044)
Vaccinia Capping Enzyme (Cat. #2460A/B)
mRNA Cap 2'-O-Methyltransferase (Cat. #2470A/B)

In-Fusion is a registered trademark of Takara Bio USA, Inc.
IVTpro and PrimeCap are trademarks of Takara Bio Inc.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc. If you require licenses for other use, please contact us by phone at +81 77 565 6972 or from our website at www.takarabio.com. Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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BspQ I, HQ *	1 μ l
10×BspQ I Buffer, HQ	10 μ l
基質 DNA	≤ 10 μ g
Nuclease-free water	up to 100 μ l

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