# 0.5-10 kb ssRNA Ladder Marker

**Code No. 3417A**      **Size: 50  $\mu$ l**  
**(for 25 lanes)**

**Supplied Reagents:**  
**RNA Loading Buffer**      **1 ml**  
**6X Dye**      **1 ml**

**Form:** 1 mM Sodium Citrate, pH 6.4

**Description:**

This marker consists of 8 single-stranded RNA (ssRNA) fragments ranging from 0.5 kb to 10 kb.

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**Storage:** -80°C

**Usage:**

Used as a molecular size marker for long ssRNA such as mRNA and lncRNA (long non-coding RNA) in gel electrophoresis.

**RNA Loading Buffer: (Can be stored at -20°C after opening)**

96% Formamide  
20 mM EDTA, pH 8.0

**6X Dye: (Can be stored at -20°C or 4°C after opening)**

0.5% Orange G  
1 mM EDTA, pH 8.0

**Note:**

1. Store at -80°C quickly after use.
2. Do not store ssRNA Ladder Marker after adding RNA Loading Buffer or 6X Dye.
3. Extra precaution should be taken during experiments in order to prevent RNase contamination. Wear clean disposable gloves, and use tubes and micropipette chips exclusively for RNA experiments. Do not perform experiments which use RNase, such as plasmid preparation in the same area.
4. Depending on the sequence or volume of sample ssRNA, the fragments may appear at slightly different lengths from this ladder marker.

**Application Example:**

1. Preparation of RNA marker solution.

0.5-10 kb ssRNA Ladder Marker	2 $\mu$ l
RNA Loading Buffer	3 $\mu$ l
6X Dye	1 $\mu$ l
<b>Total</b>	<b>6 <math>\mu</math>l</b>

2. Denature at 65°C for 5 min and immediately cool on ice.
3. Load 6  $\mu$ l\*1 of the prepared RNA marker solution on non-denaturing or denaturing agarose gel, and run electrophoresis.\*2
4. Stain the gel if necessary and photograph the gel under UV irradiation.\*3

\*1 Prepare the required amount of RNA marker solution depending on the well-size of the agarose gel.

\*2 Run until the dye front (Orange G) has migrated to the tip of the gel or run out.

\*3 There are several ways for staining RNA. When performing Ethidium bromide (EtBr) staining, there are three types of methods.

Prestain

Add EtBr to the agarose gel in advance to a final concentration of 1  $\mu$ g/ml.

Poststain

After electrophoresis, stain the gel with running buffer containing EtBr at a final concentration of 1  $\mu$ g/ml, followed by washing with the EtBr-free buffer for 5 - 10 min.


Direct stain

Add EtBr to the sample in advance to a final concentration of 12.5  $\mu$ g/ml.

[ Electrophoresis Examples ]

Left : 2% non-denaturing (TAE) agarose gel (Prestain)

Right: 1.2% denaturing (Formaldehyde/MOPS) agarose gel (Direct stain)

**Note**

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc.

If you require licenses for other use, please contact us by phone at +81 77 565 6972 or from our website at [www.takarabio.com](http://www.takarabio.com).

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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