EmeraldAmp[®] PCR Master Mix

Code No. RR300A Size: 1 ml x 4
(for 160 PCR reactions)

Supplied Reagent:
dH₂O 1 ml x 4

Description:

EmeraldAmp PCR Master Mix (2X Premix) is a premixed, 2 - fold concentrated PCR mixture composed of PCR enzyme for Hot Start PCR, reaction buffer and dNTP Mixture. In addition, this product also contains blue dye and yellow dye markers and a material for increasing specific gravity, useful in electrophoresis. Operation of PCR experiment is greatly simplified with using this product, in the points that PCR reaction can be done by only adding with primers and template DNA and that PCR product can be applied directly into electrophoresis analysis. The PCR reaction mixture shows vivid emerald green, therefore, its application into electrophoresis gel is easily performed.

This product is very useful for examining of an insert fragment by colony-directed PCR of *E. coli* recombinant cells, and also applicable to amplify GC-rich region of DNA. Furthermore, it is possible to amplify 10 kb of DNA fragment if used human genomic DNA as a template.

Storage: -20°C (at 4°C for 3 months)

(Note) Recommend to store at 4°C in case of used frequently, because the activity of the Master Mix may decrease by repeated freeze-thawing. Please use after gently mixing well, and briefly centrifuging.

Application:

- DNA amplification by PCR
- Colony PCR

Quality Control Data :

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

PCR product:

As most PCR products amplified with EmeraldAmp PCR Master Mix have one A added at 3'-termini, the obtained PCR product can be directly used for cloning into T-vector. Also it is possible to clone the product in blunt-end vectors after blunting and phosphorylation of the end.

Dye marker at electrophoresis:

If 5 μ l of the reaction mixture is used for electrophoresis with 1% Agarose L03 [TAKARA] (Cat. #5003), the blue dye marker is detected near 3 - 5 kb and the yellow is below 50 bp. Those dyes have absorptions at around 260 nm and 420 nm, respectively, so it is recommended to remove the dye by cutting out the gel or extracting DNA by NucleoSpin Gel and PCR Clean-up (Cat. #740609.10/.50/.250).

General reaction mixture for PCR (Total 50 μ l):

EmeraldAmp PCR Master Mix (2X Premix)	25 μ l				
Forward Primer	0.2 μ M (final conc.)				
Reverse Primer	0.2 μ M (final conc.)				
Template	<table border="0"><tr><td>┌ human genomic DNA</td><td>50 - 100 ng</td></tr><tr><td>└ plasmid DNA</td><td>1 pg - 1 ng</td></tr></table>	┌ human genomic DNA	50 - 100 ng	└ plasmid DNA	1 pg - 1 ng
┌ human genomic DNA	50 - 100 ng				
└ plasmid DNA	1 pg - 1 ng				
dH ₂ O	up to 50 μ l				

PCR Conditions (example):

3 step PCR (- 6 kb):	98°C	10 sec	} 30 cycles
	55°C or 60°C *1	30 sec	
	72°C	1 min/kb	
2 step PCR (6 kb -):	98°C	10 sec	} 30 cycles
	68°C	1 min/kb *2	

- * 1 Annealing temperature of 60°C may be better for amplification of GC-rich regions.
- * 2 Amplified product may become smear if the extension time is set over 15 min.

(Note) Denaturation conditions vary depending on the thermal cycler and tubes used for PCR. The recommendation is for 5 - 10 sec at 98°C, or 20 - 30 sec at 94°C.

EmeraldAmp is a registered trademark of Takara Bio Inc.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc. If you require licenses for other use, please contact us by phone at +81 77 565 6972 or from our website at www.takarabio.com. Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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