# Pyrithiamine-resistant Transformation Systems

## Transformation of *A. oryzae* through Protoplast-PEG method

Sterilize all instruments and solutions to be used, before starting the protocol.

All operations must be performed under sterile condition.

### [Preparation]

• Spore suspension :

Collect the spores of *A. oryzae* on a plate and suspend in 10 ml of 0.1% Tween 80 and 0.8% NaCl. Filter the suspension with a glass filtration (3G2) to remove the hypha and the agar debris, and recover the filtered spore solution.

Centrifuge at 3,000 rpm for 5 min to collect the spores, and discard the supernatant. Wash the spores with 10 ml of 0.1% Tween 80 twice. Suspend the spores in a proper amount of sterile purified water.

• CD medium (per 1 liter) :

NaNO <sub>3</sub>	6.0 g
KCl	0.52 g
KH <sub>2</sub> PO <sub>4</sub>	1.52 g
Glucose	10.0 g *1
1 M MgSO <sub>4</sub> · 7H <sub>2</sub> O	2 ml *2
Trace elements solution *3	1 ml
Agar	20.0 g)

Adjust to pH 6.5 with 1 N KOH.

- \* 1 MgSO<sub>4</sub> · 7H<sub>2</sub>O should be autoclaved separately and mixed later.
- \* 2 When a strain shows very bad growth, culture at the higher glucose concentration around up to 20.0 gram.
- \* 3 Trace elements solution (per 1 liter)
 

FeSO <sub>4</sub> · 7H <sub>2</sub> O	1.0 g
ZnSO <sub>4</sub> · 7H <sub>2</sub> O	8.8 g
CuSO <sub>4</sub> · 5H <sub>2</sub> O	0.4 g
Na <sub>2</sub> B <sub>4</sub> O <sub>7</sub> · 10H <sub>2</sub> O	0.1 g
(NH <sub>4</sub> ) <sub>6</sub> Mo <sub>7</sub> O <sub>24</sub> · 4H <sub>2</sub> O	0.05 g

• CD selective medium :

CD medium added with 0.8 M NaCl and 0.1 μg/ml Pyrithiamine (PT).

• CD soft agar selective medium :

CD selective medium having 0.5% agar.  
Keep it at approximately 50°C.

• Protoplast forming solution : Sterilize by filtering.

20 mg/ml Yatalase (Cat. #T017), 0.8 M NaCl,  
10 mM Na phosphate buffer (pH 6.0).

• Protoplast washing solution : 0.8 M NaCl

• Solution 1 :

0.8 M NaCl, 10 mM CaCl<sub>2</sub>, 10 mM Tris-HCl (pH 8.0)

• Solution 2 : Sterilize by filtering.

40% (w/v) PEG4000, 50 mM CaCl<sub>2</sub>, 50 mM Tris-HCl (pH 8.0)

### [Protocol]

- 1) Inoculate the spore suspension of *A. oryzae* in 100 ml of CD liquid medium. \*4 Culture by shaking at 30°C for 20 hours. \*5 (Shake vigorously in a triangle flask with baffle.)
  - 2) Filter with a glass filtration (3G1) to collect the hypha. Wash with sterile purified water and remove water from the hypha fully by pressing with a spatula.
  - 3) Suspend proper amount of hypha in protoplast forming solution contained in a 50 ml polypropylene centrifuge tube.
  - 4) Shake gently at 30°C and form protoplast. Confirm the formation of protoplast by observing the status with a microscope every 30 min. (Usually the reaction is completed within 2 hours.)
  - 5) Filter with a glass filtration (3G2). Centrifuge the filtered solution at 2,000 rpm for 5 min to collect protoplast.
  - 6) Wash the protoplast twice with 0.8 M NaCl.
  - 7) Suspend the protoplast in Solution 1 to have the concentration of  $2 \times 10^8$ /ml. Then add the 20% volume of Solution 2 and suspend gently. As Solution 2 has high viscosity, suspend the solution fully.
  - 8) Add plasmid (in less than 20 μl) into 0.2 ml of protoplast suspension contained in polypropylene centrifuge tube (larger than 10 ml). Plasmid can be added at maximum 20 μg.
  - 9) Leave to stand on ice for 30 min.
  - 10) Add 1 ml of Solution 2. Suspend gently.
  - 11) Leave to stand at room temperature for 15 min.
  - 12) Add 8.5 ml of Solution 1. Suspend gently.
  - 13) Centrifuge the solution to collect the protoplast. Discard the supernatant as much as possible. Suspend the collected protoplast in 0.2 ml of Solution 1.
  - 14) Mix the protoplast suspension in 5 ml of CD soft agar selective medium. \*4 Quickly overlay the protoplast/soft agar on CD selective plate (90 mm-diameter) \*4 to be spread uniformly.
  - 15) Culture at 30°C for 5 - 7 days.
- \* 4 Use CD medium without thiamine. Under the existence of thiamine, the selection efficiency by PT drastically decreases.
  - \* 5 Some strains may allow good results after cultured for 40 - 48 hours by shaking.

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